

StairFit

Fast, safe stair training for rehabilitation

A compact set of adjustable training stairs used by physiotherapists to rebuild a patient's confidence and strength on steps before they return home.

Key features

Adjustable step height

Step height and depth can be adjusted in small increments, allowing therapists to tailor difficulty to each patient's recovery stage.

Integrated handrails

Height-adjustable rails on both sides provide confidence and support during early-stage training sessions.

V2: foldable frame

The second-generation StairFit V2 folds flat for storage and transport between clinics, wards, or home-visit sessions.

Non-slip certified surfaces

Every tread uses a certified non-slip surface tested to clinical safety standards for rehabilitation equipment.

At a glance

Step height range	10 - 20 cm, adjustable
Number of steps	3 (configurable)
Handrail height	Adjustable, both sides
Version	V2 - foldable frame
Certification	Clinical non-slip standard

Standard vs. optional

Standard

- Adjustable-height handrails, both sides
- Certified non-slip tread surfaces
- 3-step configurable frame
- Locking wheels for safe transport

Optional

- Foldable V2 frame upgrade
- Additional step modules (4th/5th step)
- Therapist clipboard tray attachment
- Padded handrail grips

Used in

